

## CT Protocol

### A. Adult OMI Protocol

- a. Performed on decedents about age 8 and older
- b. Full body adult CT acquired in three acquisitions over three body planes
- c. Head/Neck/Upper extremity protocol
  - i. Arms crossed over abdomen
  - ii. If unable to break rigor mortis, see Protocol B.
  - iii. Scan above skull vertex to below finger tips- one acquisition
  - iv. Protocol Scan Parameters

kVp	120
mAs	300
Scan length	600-800 mm
Scan FOV	500-699 mm
Pitch	0.567
Collimation	16 x 0.75
Rotation Time	1.0 s
Matrix	512 x 512

- v. Reconstructions
  1. Upper extremity reconstructions- axial
    - a. 3 x 3 mm soft tissue = 190 images
    - b. 3 x 3 mm bone = 190 images
    - c. 1 x 0.5 mm soft tissue = 1200 images
    - d. 1 x 0.5 mm bone = 1200 images
  2. Head/neck reconstructions- axial
    - a. 3 x 3 mm soft tissue = 120 images
    - b. 3 x 3 mm bone = 120 images
    - c. 1 x 0.5 mm soft tissue = 670 images
    - d. 1 x 0.5 mm bone = 670 images
    - e. 2 x 1 brain = 50 images

- d. Torso protocol
  - i. Break rigor, lifting arms above head and out of field of view for scan
  - ii. Scan above clavicles through ischial tuberosities- one acquisition
  - iii. Protocol Scan Parameters

kVp	120
mAs	300
Scan length	600-800 mm
Scan FOV	350-699 mm
Pitch	0.817
Collimation	16 x 0.75
Rotation Time	1.0 s
Matrix	512 x 512

- iv. Reconstructions
  1. Torso reconstructions- axial
    - a. 3 x 3 mm soft tissue = 270 images
    - b. 3 x 3 mm bone = 270 images
    - c. 1 x 0.5 mm soft tissue = 1600 images

- d. 1 x 0.5 mm bone = 1600 images
- e. 1 x 0.5 mm Lung = 600 images

e. Lower extremities protocol

- i. Overlap scan at ischial tuberosities and scan through toes- one acquisition
- ii. Attempt to separate lower extremities during scan to avoid superimposition
- iii. Protocol Scan Parameters

kVp	120
mAs	200
Scan length	800-1000 mm
Scan FOV	500 mm
Pitch	0.942
Collimation	16 x 0.75
Rotation Time	1.0 s
Matrix	512 x 512

iv. Reconstructions

- 1. Lower extremities reconstructions- axial
  - a. 3 x 3 mm soft tissue = 320 images
  - b. 3 x 3 mm bone = 320 images
  - c. 1 x 0.5 mm soft tissue = 1900 images
  - d. 1 x 0.5 mm bone = 1900 images

**B. Decomposed Adult OMI Protocol**

- a. Performed on decedents about age 8 and older
- b. Full body adult CT acquired in three acquisitions over three body planes
- c. Performed with decedent sealed in body bag – arms remain at sides of body
- d. Head/Neck
  - i. Scan above skull vertex to carina
  - ii. Protocol Scan Parameters

kVp	120
mAs	300
Scan length	300 mm
Scan FOV	290 mm
Pitch	0.567
Collimation	16 x 0.75
Rotation Time	1.0 s
Matrix	512 x 512

iii. Reconstructions

- 1. Head/neck reconstructions- axial
  - a. 3 x 3 mm soft tissue = 120 images
  - b. 3 x 3 mm bone = 120 images
  - c. 1 x 0.5 mm soft tissue = 670 images
  - d. 1 x 0.5 mm bone = 670 images
  - e. 2 x 1 brain = 50 images

e. Torso protocol

- i. Scan above shoulders through ischial tuberosities- one acquisition

- ii. Ensure upper extremities completely imaged
- iii. Protocol Scan Parameters

kVp	120
mAs	300
Scan length	800-1000 mm
Scan FOV	699 mm
Pitch	0.817
Collimation	16 x 0.75
Rotation Time	1.0 s
Matrix	512 x 512

iv. Reconstructions

- 1. Upper extremity/Torso reconstructions- axial
  - a. 3 x 3 mm soft tissue = 300 images
  - b. 3 x 3 mm bone = 300 images
  - c. 1 x 0.5 mm soft tissue = 1800 images
  - d. 1 x 0.5 mm bone = 1800 images
  - e. 1 x 0.5 mm Lung = 600 images

f. Lower extremities protocol

- i. Overlap scan at ischial tuberosities and scan through toes- one acquisition
- ii. Attempt to separate lower extremities during scan to avoid superimposition
- iii. Protocol Scan Parameters

kVp	120
mAs	200
Scan length	800-1000 mm
Scan FOV	500 mm
Pitch	0.942
Collimation	16 x 0.75
Rotation Time	1.0 s
Matrix	512 x 512

iv. Reconstructions

- 1. Lower extremities reconstructions- axial
  - a. 3 x 3 mm soft tissue = 320 images
  - b. 3 x 3 mm bone = 320 images
  - c. 1 x 0.5 mm soft tissue = 1900 images
  - d. 1 x 0.5 mm bone = 1900 images

**C. Pediatric OMI Protocol**

- a. Performed in one acquisition if decedent is less than 36 inches long
  - i. Scan above skull vertex through toes- one acquisition
  - ii. Ensure upper and lower extremities completely images
  - iii. Protocol Scan Parameters

kVp	120
mAs	200
Scan length	400 mm
Scan FOV	180 mm
Pitch	0.938

Collimation	16 x 0.75
Rotation Time	0.50 s
Matrix	512 x 512

iv. Reconstructions

1. Whole body reconstructions- axial

- a. 3 x 3 mm soft tissue = 150 images
- b. 3 x 3 mm bone = 150 images
- c. 1 x 0.5 mm soft tissue = 1200 images
- d. 1 x 0.5 mm bone = 1200 images
- e. 1 x 0.5 mm Lung = 250 images
- f. 2 x 0.5 mm brain = 40 images

b. Performed in two acquisitions if decedent is over 36 inches long

i. Upper extremities

1. Scan above skull vertex through ischial tuberosities- one acquisition
2. Ensure upper extremities completely imaged
3. Protocol Scan Parameters

kVp	120
mAs	200
Scan length	600-800 mm
Scan FOV	180 mm
Pitch	0.938
Collimation	16 x 0.75
Rotation Time	0.50 s
Matrix	512 x 512

4. Reconstructions

a. Skull through ischial tuberosities reconstruction- axial

- i. 3 x 3 mm soft tissue = 150 images
- ii. 3 x 3 mm bone = 150 images
- iii. 1 x 0.5 mm soft tissue = 1600 images
- iv. 1 x 0.5 mm bone = 1600 images
- v. 1 x 0.5 mm Lung = 350 images
- vi. 2 x 0.5 mm brain = 40 images

ii. Lower extremities

1. Overlap scan at ischial tuberosities and scan through toes- one acquisition
2. Attempt to separate lower extremities during scan to avoid superimposition
3. Protocol Scan Parameters

kVp	120
mAs	125
Scan length	400-600 mm
Scan FOV	300 mm
Pitch	0.942
Collimation	16 x 0.75
Rotation Time	1.0 s
Matrix	512 x 512

4. Reconstructions

a. Lower extremities reconstructions- axial

- i. 3 x 3 mm soft tissue = 125 images
- ii. 3 x 3 mm bone = 125 images
- iii. 1 x 0.5 mm soft tissue = 1200 images
- iv. 1 x 0.5 mm bone = 1200 images

**D. Partial OMI Protocol**

- a. One of the OMI protocols in part
- b. The complete body is not captured